

NMR Support Site アカウント作成ガイド

※Deltaのダウンロード等は、アカウントが承認されるまでできません。
アカウント作成から承認までにひと晩程度を要するので、ご注意ください


NMR Support Siteでアカウントを作成すると、

- ・ Deltaソフトウェアのダウンロード
 - ・ ライセンスキーの発行
 - ・ 掲示板の閲覧
- などができます。

※アカウントの作成には、有効なメールアドレスが必要です。

それでは、アカウントの作成を始めましょう。

まずは「[2. Open a User Account](#)」をクリックしてください。



<https://nmrsupport.jeol.com/>

A screenshot of the JEOL Delta / NMR Support website. The page features a header with the JEOL logo and navigation links. The main content area includes a welcome message, a list of steps to get started, and various resource links. A blue arrow points from the text '2. Open a User Account' in the main text to the corresponding link on the website. The website footer contains contact information and copyright details.

JEOL
Solutions for Innovation

My Account [Social Media Icons]

Home Docs FAQs Videos

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

1. [Login to Your Existing Account](#)
You will have full content and features available for your account type.
2. [Open a User Account](#)
Registration for this site is moderated - to gain full access, be sure to **also submit your spectrometer information** so that we can assign you the proper authorization for software downloads and/or license key generation.

Once your account registration and spectrometer information is reviewed, you will be granted the appropriate access to other areas of the site.

Read Me

Title

- [NMR Support Read Me \(English\)](#)
- [NMR Support Read Me 日本語 \(Japanese\)](#)

Other Resources

- [JEOL USA NMR](#)
Products, systems and applications that illustrate JEOL's NMR Spectroscopy expertise.
[Learn More >](#)
- [JEOL USA NMR RESOURCES](#)
Documentation in support of your JEOL product.
[Learn More >](#)
- [JEOL USA WEBINARS](#)
View all available webinars.
[Learn More >](#)
- [JEOL USA YOUTUBE CHANNEL](#)
Product overview presentations, deep-dive tutorials, and trade show demonstrations.
[Learn More >](#)

If you have comments or questions regarding this web site, please contact the **site administrator**.
DO NOT send questions about Delta operation or installation. DO NOT request Delta License Keys from the administrator.

JEOL USA, Inc.
11 Dearborn Road
Peabody, MA 01960

Please direct questions or comments about Delta or JEOL NMR spectrometers to the discussion boards.

© Copyright 2007-2020 by JEOL USA, Inc. [Terms of Use](#) | [Privacy Policy](#)

The image shows a web browser window displaying the registration page for a JEOL account. The page title is "Register for a JEOL Delta/NMR Su" and the URL is "https://nmrsupport.jeol.com/Home/Account-Registration". The JEOL logo is in the top left, and social media icons for Facebook, Instagram, LinkedIn, Twitter, and YouTube are in the top right. The main content is a "Sign up" form with the following fields: a name field containing "Katsuo Asakura", an email field containing "katsuo.asakura.jeol@gmail.com", a password field with masked characters ".....", a checkbox for "I accept the Terms of Use & Privacy Policy" which is checked, and a green "Sign up" button. Below the button is a link: "Already have an account? [Sign in](#)". Red arrows point from the Japanese text on the right to each of these form elements.

名前を設定します
(アルファベットで設定することを推奨します)

メールアドレスを入力します

好きなパスワードを設定します
(後から変更することも可能です)

チェックします
(利用許諾条件の承認です)

上記を入力したらクリックしましょう

If you have comments or questions regarding this web site, please contact the **site administrator**.

DO NOT send questions about Delta operation or installation. DO NOT request Delta License Keys from the administrator.

JEOL USA, Inc.
11 Dearborn Road
Peabody, MA 01960

Please direct questions or comments about Delta or JEOL NMR spectrometers to the discussion boards.



ユーザー情報の入力画面が現れます

The screenshot shows a web browser window with the URL <https://nmrsupportjeol.com/Home/Account-Registration>. The page features the JEOL logo (Solutions for Innovation) in the top left. A user profile card for "Katsuo Asakura" is displayed, with a blue "Profile" button. A yellow warning banner at the top of the profile card reads "Please update your profile before continuing." Below this, a "Profile - required" dropdown menu is open, showing a list of profile options. The form fields include:

- First (Given) Name:** [Text input field]
- Family (Last) Name:** [Text input field]
- Job Title:** [Dropdown menu, currently set to "< Not Specified >"]
- Company or University Name: (No Abbreviations Please!)** [Text input field]
- Address:** [Text input field]
- Address (2):** [Text input field]
- Address (3):** [Text input field]
- City:** [Text input field]

An "Update" button is located at the bottom right of the profile card. At the bottom of the page, there is a footer with the text: "If you have comments or questions regarding this web site, please contact the site administrator."

Register for a JEOL Delta/NMR Support Account

https://nmrsupport.jeol.com/Home/Account-Registration

JEOL Solutions for Researchers

Katsuo Asakura

Profile

Please update your profile before continuing

Profile - required

First (Given) Name: Katsuo

Family (Last) Name: Asakura

Job Title: Faculty

Company or University Name: (No Abbreviations Please) Japan Electron Optics Laboratory University

Address: 3-1-2 Musashi-no

Address (2):

Address (3):

City: Akishima

Country: Japan (JP)

Region: Tōkyō

Postal Code: 1968558

Telephone: 042-543-1111

How Did You Find Out about Delta? Friend

Profile - optional

Preferences

Spectrometer Information

Update

名前

苗字

職業など

所属機関（大学名等）

アルファベット頭文字の略称は使用しないでください

市町村以下の住所

市町村名

国（Japan (JP)を選択します）

都道府県を選択

郵便番号

電話番号

Deltaを知った経緯を選択

上記をすべて入力したらクリックしましょう

装置情報の入力画面が現れます

Register Spectrometer

https://nmrsupportjeol.com/Home/Spectrometer-Registration

Katsuo Asakura

Home Docs Forums FAQs Videos

Spectrometer and System Information, Page 2 of 2

Your registration is NOT complete until you supply the information requested below.

Spectrometer Brand

Spectrometer Field

Workstation

Operating System

Operating System Version

Do you need Delta for:

Data Processing Only

Data Processing and JEOL NMR Machine Control (JEOL NMR Spectrometers Only)

Are you a JEOL employee or representative? No Yes

Register Spectrometer

https://nmrsupport.jeol.com/Home/Spectrometer-Registration

Katsuo Asakura

Home Docs Forums FAQs Videos

Spectrometer and System Information, Page 2 of 2

Your registration is NOT complete until you supply the information requested below.

Spectrometer Brand: JEOL

JEOL Spectrometer Type: ECZLS

Spectrometer Field: 400

Workstation: Intel PC

Operating System: Windows

Operating System Version: 11

Do you need Delta for:

Data Processing Only

Data Processing and JEOL NMR Machine Control (JEOL NMR Spectrometers Only)

Are you a JEOL employee or representative? No Yes

Submit

普段使っているNMR装置のメーカー

普段使っているNMR装置のモデル名
(メーカー名を選ぶと現れます)

普段使っているNMR装置の磁場強度

普段使っているPCの種類

普段使っているPCのOS

普段使っているPCのOSのバージョン

Deltaの使用目的

※必ず「Data Processing Only」を選びます


JEOL社員かどうか

※必ず「No」を選びます

上記をすべて入力したらクリックしましょう


メールの受信について注意喚起する画面が現れます

この画面では、何もクリックせずに、メールが届いたかどうか確認しましょう。
見当たらないときは迷惑メールフォルダも確認しましょう。



しばらく待っても**メール**
が届かない場合のみ、
ここをクリックして、
メールアドレスを修正し
ます。


こんなメールが届きます。 (アカウントの作成を開始した段階で届いています)



**※この時点ではまだ承認されていないので、
ダウンロード等ができません**

こんなメールが届きます。

(装置情報を入力・送信した段階で届きます)



The screenshot shows a Gmail interface with an email from "JEOL Webmaster" (nmr-admin@jeol.com) dated March 12, 2026, at 3:01 AM. The email subject is "Spectrometer Information submitted from JEOL USA NMR Support for Asakura, Katsuo(katsuo.asakura.jeol@gmail.com)". The body of the email contains the following information:

Thursday, March 12, 2026 3:01 AM
Dear Katsuo Asakura:
Thank you for providing the following spectrometer information.

Username: katsuo.asakura.jeol@gmail.com
Job Title: Faculty
Institution or Company: Japan Electron Optics Laboratory University
Email: katsuo.asakura.jeol@gmail.com

Spectrometer Brand: JEOL
Spectrometer Type: ECZS
Spectrometer Field: 400

Workstation: Intel PC
Operating system: Windows
Operating system version: 11


Delta Type: Data Processing Only
Spectrometer S/N:

JEOL Affiliated: No
Company / Distributor Name:

If appropriate, you will be contacted after your submission is reviewed.
Thank you for your interest in Delta.
JEOL NMR Staff

At the bottom of the email, there are buttons for "返信" (Reply), "転送" (Forward), and a smiley face icon. A large red text overlay is present at the bottom of the screenshot, stating: "※この時点ではまだ承認されていないので、ダウンロード等ができません" (※At this point, since it has not been approved yet, downloading, etc., is not possible).

一晩くらいでこんなメールが届きます。 (アカウントが承認された段階で届きます)








The screenshot shows a Gmail interface with a message titled "JEOL Delta / NMR Support Site User Account Update". The sender is "nmr-admin@jeol.com". The email content reads: "Dear Katsuo Asakura, Your user account at JEOL Delta / NMR Support Site has been recently updated to include access to the following: Area: DELTA Description: Data Processing Only Account Thank you, we appreciate your support. JEOL Delta / NMR Support Site". Below the email content are buttons for "返信" (Reply), "転送" (Forward), and a "受信トレイ" (Inbox) icon. The left sidebar shows navigation options like "作成" (Compose), "受信トレイ" (Inbox), "スター付き" (Starred), "スヌーズ中" (On mute), "送信済み" (Sent), "下書き" (Drafts), "購入" (Purchased), and "もっと見る" (See more). The bottom left corner has an "アップグレード" (Upgrade) button.

※無事に承認されたので、ダウンロード等ができるようになりました！


NMR Support Siteにサインインしましょう！

JEOL Delta / NMR Support Site

https://nmrsupport.jeol.com

My Account     

Home Docs FAQs Videos



Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1 Login to Your Existing Account**
You will then see content and features available for your account type.
- 2 Open a User Account**
Registration to this site is moderated - to gain full access, be sure to **also submit your spectrometer information** so that we can assign you the proper authorization for software downloads and for license key generation.





Once your account registration and spectrometer information is reviewed, you will be granted the appropriate access to other areas of the site.

Read Me

Title

- [NMR Support Read Me \(English\)](#)
- [NMR Support Read Me 日本語 \(Japanese\)](#)






Other Resources

 JEOL USA NMR Products, systems and applications that illustrate JEOL's NMR Spectroscopy expertise. Learn More >	 JEOL USA NMR RESOURCES Documentation in support of your JEOL product. Learn More >	 JEOL USA WEBINARS View all available webinars. Learn More >	 JEOL USA YOUTUBE CHANNEL Product overview presentations, deep-dive tutorials, and trade show demonstrations. Learn More >
--	---	--	--

こんどは「**1. Login to Your Existing Account**」をクリックしてください。

JEOL Delta / NMR Support Site

https://nmrsupportjeol.com/Home/SignIn?returnurl=%2f

My Account     

Home Docs FAQs Videos

Sign in

Keep me signed in


[Lost your password?](#)

[Don't have an account? Sign up](#)

If you have comments or questions regarding this web site, please contact the **site administrator**.
DO NOT send questions about Delta operation or installation. DO NOT request Delta License Keys from the administrator.

JEOL USA, Inc.
11 Dearborn Road
Peabody, MA 01960

Please direct questions or comments about Delta or JEOL NMR spectrometers to the discussion boards.



The screenshot shows the JEOL Delta / NMR Support website. The top navigation bar includes links for Home, Software, License Keys, Docs, Forums, FAQs, and Videos. The main content area features a large image of a hand holding a pipette tip. Below the image, the heading "Welcome to JEOL Delta / NMR Support" is followed by a paragraph explaining the site's purpose. A section titled "To get started, simply:" lists two steps: "1. Login to Your Existing Account" and "2. Open a User Account". The "Open a User Account" step includes a note that registration is moderated and requires submitting spectrometer information. Below this, a paragraph states that once registration and spectrometer information is reviewed, appropriate access will be granted. A red rounded rectangle highlights the "Open a User Account" step and the subsequent paragraph. Below the highlighted area, the same heading and paragraph are repeated, followed by a list of seven options: "1. Download Delta Software", "2. Obtain a License Key", "3. Document Downloads", "4. Discussion Forum", "5. FAQs", "6. Manage / Change User Profile", and "7. Change / Reset Password". The "Read Me" section offers links for English and Japanese versions. The "Other Resources" section includes links to JEOL USA NMR products, NMR resources, webinars, and a YouTube channel.

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1. Login to Your Existing Account**
You will then see content and features available for your account type.
- 2. Open a User Account**
Registration to this site is moderated - to gain full access, be sure to **also submit your spectrometer information** so that we can assign you the proper authorization for software downloads and/or license key generation.

Once your account registration and spectrometer information is reviewed, you will be granted the appropriate access to other areas of the site.

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1. Download Delta Software**
Depending on your account type, Delta is available for data processing only, or for full JEOL NMR spectrometer control. The current and older versions of Delta are available for download.
- 2. Obtain a License Key**
Submit the Delta Machine ID of your workstation to obtain a matching Delta License Key to enable the software.
- 3. Document Downloads**
Download Release Notes, Delta Tutorials, User Guides, and Other Documents
- 4. Discussion Forum**
Discussion Forums Related to Delta and JEOL NMR Spectrometers
- 5. FAQs**
Frequently Asked Questions related to Delta and JEOL NMR Spectrometers.
- 6. Manage / Change User Profile**
Manage contact and spectrometer information.
- 7. Change / Reset Password**
Change your Password, Click "Send Reset Link"

Read Me

Title
NMR Support Read Me (English)
NMR Support Read Me 日本語 (Japanese)

Other Resources

- JEOL USA NMR**
Products, systems and applications that illustrate JEOL's NMR Spectroscopy expertise.
[Learn More >](#)
- JEOL USA NMR RESOURCES**
Documentation in support of your JEOL product.
[Learn More >](#)
- JEOL USA WEBINARS**
View all available webinars.
[Learn More >](#)
- JEOL USA YOUTUBE CHANNEL**
Product overview presentations, deep-dive tutorials, and trade show demonstrations.
[Learn More >](#)

この部分が表示されていれば
アカウントの作成は完了です

Deltaのダウンロード

Deltaソフトウェアは、ダウンロードしたインストーラでPCにインストールします。インストーラファイルは、複製したり、他のPCに移動したりしても構いません。複数台のPCにインストールすることも問題ありません。

※Deltaを使うときにインターネットに接続している必要はありません。

それでは、Deltaのダウンロードを始めましょう。

「[1. Download Delta Software](#)」をクリックしてください。

JEOL Delta / NMR Support Site

Home Software License Keys Docs Forums FAQs Videos

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1. Login to Your Existing Account**
You will then see content and features available for your account type.
- 2. Open a User Account**
Registration to this site is moderated -- to gain full access, be sure to **also submit your spectrometer information** so that we can assign you the proper authorization for software downloads and/or license key generation.

Once your account registration and spectrometer information is reviewed, you will be granted the appropriate access to other areas of the site.

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1. Download Delta Software**
Depending on your spectrometer type, Delta is available for data processing only, or for full JEOL NMR spectrometer control. The current and older versions of Delta are available for download.
- 2. Obtain a License Key**
Submit the Delta Machine ID of your workstation to obtain a matching Delta License Key to enable the software.
- 3. Document Downloads**
Download Release Notes, Delta Tutorials, User Guides, and Other Documents
- 4. Discussion Forum**
Discussion Forums Related to Delta and JEOL NMR Spectrometers
- 5. FAQs**
Frequently Asked Questions related to Delta and JEOL NMR Spectrometers.
- 6. Manage / Change User Profile**
Manage contact and spectrometer information.
- 7. Change / Reset Password**
Change your Password. Click "Send Reset Link"

Read Me

Title

NMR Support Read Me (English)

NMR Support Read Me 日本語 (Japanese)

Other Resources

- JEOL USA NMR
Products, systems and applications that illustrate JEOL's NMR Spectroscopy expertise.
Learn More >
- JEOL USA NMR RESOURCES
Documentation in support of your JEOL product.
Learn More >
- JEOL USA WEBINARS
View all available webinars.
Learn More >
- JEOL USA YOUTUBE CHANNEL
Product overview presentations, deep-dive tutorials, and trade show demonstrations.
Learn More >

Software

https://nmrsupportjeol.com/Software

Katsuo Asakura

Home Software License Keys Docs Forums FAQs Videos

Delta Software Download

CAREFULLY READ THE FOLLOWING TERMS AND CONDITIONS. YOU AGREE TO ALL TERMS AND CONDITIONS BY INSTALLING THE DELTA NMR SOFTWARE PACKAGE OR USING ANY PART THERE OF.

License: With the payment of any and all applicable License Fees upon the delivery of Delta NMR Software ("Software"), JEOL USA, INC. ("Licensor") grants to the Licensee a License to use without exclusively a single non-transferable, non-sublicensable copy of the Software on a single CPU. No title, ownership of, or rights to the software is hereby transferred to the licensee. The Software and its ancillary documentation may not be modified, reverse engineered, sublicensed, distributed, sold as a service on a time-sharing system, or used for any use other than the licensee's own internal use. This Software is composed of media, computer programs, and supporting documentation. In all cases the Software and its associated documentation may not be copied in whole or in part by any means or in any way disclosed to a third party. Upon License termination all Media, Software and Documentation must be removed from the CPU's and premises of the Licensee and returned to the licensor. The obligations of non-replication and non-disclosure set forth above shall survive any termination of this license.

Warranty (License Fee Based): The Licensor warrants that the software will process the supplied test data in accordance with the ancillary documentation within the limitations of the license agreement and those specified in this document for a period of one year from the date of payment of any applicable fees. The Licensor does not assume any responsibility for the use or reliability of its software or for any losses or damages which may result from its use due to negligence or otherwise, even if the Licensor has been advised, known, or should have known, of the possibility of such damages. The Licensor may or may not, during the limited warranty period, at its option agree to rectify, replace, refund the purchase price, or terminate the license agreement invoking the termination clauses listed above at its option as a response to any defects in workmanship, omissions in documentation, or oversights in design not caused or contributed to by operation of the software contrary to the terms of the operating manual and/or the term of license.

Warranty (No License Fee Based): If there are no applicable fees the Licensor makes NO warranty as to selection of this product to achieve any intended purposes, for the proper installation and/or use of this product or for verifying the results obtained from use of this product.

General: The License states the complete agreement of the parties with respect to the rights and privileges associated with the Software. This agreement supersedes all prior Licenses and agreements both written and oral. If any provision in this license is determined to be illegal or unenforceable, this License shall remain valid and all other provisions shall be enforceable in accordance with their respective terms, as fully as though such illegal or unenforceable provisions had never been a part of this License. Failure on the part of the Licensor to enforce any provision in this license in any way shall not be deemed to be a waiver of such provision or any other provision. The laws of the state of Delaware shall govern all terms and conditions of this license.

Your Environment:

Operating System

- Windows
- MAC OS
- Linux
- IRIX
- AIX
- Solaris

Workstation

- Intel PC
- MAC PPC
- MAC Intel
- MAC Apple Silicon
- SPARC
- IBM RS/6000
- Sun

Package Type

- Data Processing Only
- Full Package with Instrument Control

Windows version

11

Accept Terms and Conditions

If you have comments or questions regarding this web site, please contact the **site administrator**.
DO NOT send questions about Delta operation or installation. DO NOT request Delta License Keys from the administrator.

JEOL USA, Inc.
11 Dearborn Road
Peabody, MA 01960

Please direct questions or comments about Delta or JEOL NMR spectrometers to the discussion boards.

© Copyright 2007-2026 by JEOL USA, Inc. Terms of Use | Privacy Policy

PCのOSの種類

OSのバージョン

PC (CPU) の種類

ライセンスの種類

※必ずData Processing Onlyを選択します

上記を入力したらクリックしましょう

JEOL Delta / NMR Support Site > x +

https://nmrsupportjeol.com/Software/Delta-Software-Downloads-for-Windows

Katsuo Asakura | f | @ | in | t | v

Home Software License Keys Docs Forums FAQs Videos

Delta NMR Software Version Compatibility

All NMR data opened or processed with Delta V5.x CANNOT be opened with Delta V4.3.x. See the Delta V5.x release notes for more information. JEOL is not responsible for lost or damaged NMR data when using any version of Delta.

Delta NMR data files processed and saved with Delta V5.x cannot be opened with Delta V4.3.6.

The Delta V4.3.x and V5.0.1 restriction related to opening data had been removed in Delta V5.0.1 and later versions of Delta. The restriction for processed data has not changed.

Delta Software Downloads for Windows

Package Type	Workstation	Package Name	Operating System	OS Versions	Delta Version	Delta NMR software file
Data Processing Only Full Package with Instrument Control	Intel PC	Delta for Windows v6.5.0	Windows	Windows 10 and newer	v6.5.0	Delta_V650_Windows.zip
Data Processing Only Full Package with Instrument Control	Intel PC	Delta for Windows v6.4.0	Windows	Windows 10 64 Bit and Windows11	v6.4.0	Delta_V640_Windows.zip
Data Processing Only Full Package with Instrument Control	Intel PC	Delta for Windows v6.3.0	Windows	Windows 10 64 Bit and Windows11	v6.3.0	Delta_v630_Windows.zip
Data Processing Only Full Package with Instrument Control	Intel PC	Delta for Windows V6.2.0	Windows	Windows 10 64 Bit and Windows11	V6.2.0	Delta_v620_Windows.zip
Data Processing Only Full Package with Instrument Control	Intel PC	Delta for Windows V6.1.0	Windows	Windows 10 64 Bit and Windows11	V6.1.0	Delta_v610_Windows.zip
Data Processing Only	Intel PC	Delta for Windows V6.0.0	Windows	Windows 10 64 Bit	V6.0.0	delta-6.0.0-windows-x64-installer-release.exe
Data Processing Only	Intel PC	Delta for Windows V5.3.3	Windows	Windows 10 32 and 64 Bit	V5.3.3	Delta_V5.3.3_Windows.zip
Data Processing Only	Intel PC	Delta for Windows v5.3.1	Windows	Current	5.3.1	Delta_V5.3.1_Windows.zip
Data Processing Only	Intel PC	Delta for Windows v5.3.0	Windows	Current	5.3.0	Delta_V5.3.0_Windows.zip
Data Processing Only	Intel PC	Delta for Windows v5.2.1	Windows	Current	5.2.1	Delta_V5.2.1_Windows.zip

1 2 3 >>


インストーラファイル

ダウンロードしたインストーラファイルはZIP圧縮されていますので、解凍してインストールしてください。


ライセンスキーの発行

インストールしたDeltaは15分間しか使えません
無料のライセンスキーを設定して、制限を解除しましょう


インストールしたDeltaを起動すると、
マスターコンソールと呼ばれるウィンドウが開きます。



オプションメニューの「ライセンスキー…」を選択します



ライセンスキーの設定ウィンドウが開きます。



Machine IDが表示されていますので、
NMR Support Siteで入力し、ライセンスキーを発行します
※Machine IDはコンピューター1台ごとに異なります

それでは、ライセンスキーの発行を始めましょう。

「 2. Obtain a License Key 」をクリックしてください。

JEOL Delta / NMR Support Site

https://nmsupport.jeol.com

Katsuo Asakura

Home Software License Keys Docs Forums FAQs Videos

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1. Login to Your Existing Account**
You will then see content and features available for your account type.
- 2. Obtain a License Key**
Registration to this site is moderated -- to gain full access, be sure to **also submit your spectrometer information** so that we can assign you the proper authorization for software downloads and/or license key generation.

Once your account registration and spectrometer information is reviewed, you will be granted the appropriate access to other areas of the site.

Welcome to JEOL Delta / NMR Support

The JEOL Delta / NMR support web site provides support for users interested in the Delta NMR Software package and JEOL's NMR Spectrometers. The web site includes Delta NMR software kits for NMR data processing only and kits for JEOL NMR Spectrometer control with data processing. Areas for Delta License Key generation, discussion forums, FAQ's, and documentation downloads are also available on this web site.

To get started, simply:

- 1. Download Delta Software**
Depending on your account type, Delta is available for data processing only, or for full JEOL NMR spectrometer control. The current and older versions of Delta are available for download.
- 2. Obtain a License Key**
Enter the ID of your workstation to obtain a matching Delta License Key to enable the software.
- 3. Document Downloads**
Download Release Notes, Delta Tutorials, User Guides, and Other Documents
- 4. Discussion Forum**
Discussion Forums Related to Delta and JEOL NMR Spectrometers
- 5. FAQs**
Frequently Asked Questions related to Delta and JEOL NMR Spectrometers.
- 6. Manage / Change User Profile**
Manage contact and spectrometer information.
- 7. Change / Reset Password**
Change your Password. Click "Send Reset Link"

Read Me

Title

NMR Support Read Me (English)

NMR Support Read Me 日本語 (Japanese)

Other Resources

- JEOL USA NMR**
Products, systems and applications that illustrate JEOL's NMR Spectroscopy expertise.
[Learn More >](#)
- JEOL USA NMR RESOURCES**
Documentation in support of your JEOL product.
[Learn More >](#)
- JEOL USA WEBINARS**
View all available webinars.
[Learn More >](#)
- JEOL USA YOUTUBE CHANNEL**
Product overview presentations, deep-dive tutorials, and trade show demonstrations.
[Learn More >](#)

License Keys

https://nmrsupportjeol.com/License-Keys

Katsuo Asakura

JEOL Solutions for Innovation

Home Software License Keys Docs Forums FAQs Videos

Delta License Key Generation

Terms and Conditions

CAREFULLY READ THE FOLLOWING TERMS AND CONDITIONS. YOU AGREE TO ALL TERMS AND CONDITIONS BY INSTALLING THE DELTA NMR SOFTWARE PACKAGE OR USING ANY PART THERE OF.

License: With the payment of any and all applicable License Fees upon the delivery of Delta NMR Software (Software), JEOL USA, INC. (Licensor) grants to the Licensee a License to use without exclusivity a single non-transferable, non-sublicensable copy of the Software on a single CPU. No title, ownership of, or rights to the software is hereby transferred to the licensee. The Software and its ancillary documentation may not be modified, reverse engineered, sublicensed, distributed, sold as a service on a time-sharing system, or used for any use other than the licensee's own internal use. This Software is composed of media, computer programs, and supporting documentation. In all cases the Software and its associated documentation may not be copied in whole or in part by any means or in any way disclosed to a third party. Upon License termination all Media, Software and Documentation must be removed from the CPU's and

Accept Terms Reject Terms

PackageType

Data Processing Only Full Key with Instrument Control

Workstation

Operating System

Operating System Version

Delta Machine ID

NOTE: The Delta Machine ID is generated by Delta not the operating system. It is not the MAC address.

Read Me

Title

[NMR Support Read Me \(English\)](#)

[NMR Support Read Me 日本語 \(Japanese\)](#)

If you have comments or questions regarding this web site, please contact the [site administrator](#).

JEOL USA, Inc.
11 Dearborn Road
Peabody, MA 01960

Please direct questions or comments about Delta or JEOL NMR spectrometers to the [discussion boards](#).

ソフトウェア利用許諾条件の承諾

※Accept Termsをチェックします

ライセンスの種類

※必ずData Processing Onlyを選択します

PC (CPU) の種類

PCのOSの種類

OSのバージョン

Machine ID

※Deltaで確認したIDを入力します

0 (ゼロ) と O (オー)、1 (イチ) と I (アイ) などを間違えないように気をつけましょう

上記を入力したらクリックしましょう

License Keys

https://nmrsupportjeol.com/License-Keys

Katsuo Asakura

JEOL Solutions for Innovation

Home Software License Keys Docs Forums FAQs Videos

Delta License Key Generation

Terms and Conditions

CAREFULLY READ THE FOLLOWING TERMS AND CONDITIONS. YOU AGREE TO ALL TERMS AND CONDITIONS BY INSTALLING THE DELTA NMR SOFTWARE PACKAGE OR USING ANY PART THERE OF.

License: With the payment of any and all applicable License Fees upon the delivery of Delta NMR Software (Software), JEOL USA, INC. (Licensor) grants to the Licensee a License to use without exclusivity a single non-transferable, non-sublicensable copy of the Software on a single CPU. No title to, ownership of, or rights to the software is hereby transferred to the licensee. The Software and its ancillary documentation may not be modified, reverse engineered, sublicensed, distributed, sold as a service on a time-sharing system, or used for any use other than the licensee's own internal use. This Software is composed of media, computer programs, and supporting documentation. In all cases the Software and its associated documentation may not be copied in whole or in part by any means or in any way disclosed to a third party. Upon License termination all Media, Software and Documentation must be removed from the CPU's and

Accept Terms Reject Terms

PackageType Data Processing Only Full Key with Instr

Workstation Intel PC

Operating System Windows

Operating System Version 11

Delta Machine ID VL - JN - JT -

NOTE: The Delta Machine ID is not the operating system ID.

Your Delta Machine ID is VL-JN-JT-XK-AR

Your Delta License Key is HX-QV-UV-RV-GC

This information is also being E-mailed to your address on record.

Submit

Read Me

Title

[NMR Support Read Me \(English\)](#)

[NMR Support Read Me 日本語 \(Japanese\)](#)

If you have comments or questions regarding this web site, please contact the **site administrator**. DO NOT send questions about Delta operation or installation. DO NOT request Delta License Keys from the administrator.

JEOL USA, Inc.
11 Dearborn Road
Peabody, MA 01960

Please direct questions or comments about Delta or JEOL NMR spectrometers to the discussion boards.

© Copyright 2007-2026 by JEOL USA, Inc. Terms of Use | Privacy Policy

即座にライセンスキーが発行されます

Your Delta Machine ID is
VL-JN-JT-XK-AR

こちらはMachine ID

Your Delta License Key is
HX-QV-UV-RV-GC


発行されたライセンスキー

This information is also being E-mailed to your address on record.


This information is also being E-mailed to your address on record.

・ライセンスキーはメールでも届いています

ライセンスキーを入力し、適用ボタンをクリックします



0 (ゼロ) と O (オー)、1 (イチ) と I (アイ) などを間違えないように気をつけましょう



この画面が出たら、青いボタンをクリックして設定完了です。

(この画面は表示されない場合もあります)